

**FastDigest Tai I****Thermo Scientific**

5'... A C G T#...3' 3'...#T G C A ...5' L'enzima di restrizione Tail riconosce i siti ACGT[^] e taglia a 65°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: HpyCH4IV, MaeII).



Code	Description	Packaging
FL5FD1144	FastDigest Tai I	50 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

