



FastDigest TaqI

Thermo Scientific

5'...T#C G A...3' 3'...A G C#T...5' L'enzima di restrizione TaqI riconosce i siti T[^]CGA e taglia a 65°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD0674	FastDigest TaqI	400 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

