

**FastDigest Tatl****Thermo Scientific**

5'...W#G T A C W...3' 3'...W C A T G#W...5' L'enzima di restrizione Tatl riconosce i siti W[^]GTACW e taglia a 65°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD1294	Fast Digest TatI	20 ul



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

