

**Uracil-DNA Glycosylase (1 U/ μ L)****Thermo Scientific**

Uracil-DNA Glicosilasi (UDG, UNG) catalizza l'idrolisi dei legami N-glicosilici tra uracile e zucchero nel DNA e previene il carry-over del DNA nelle reazioni di PCR.



Code	Description	Packaging
FL5EN0361	Uracil-DNA Glycosylase (1 U/ μ L)	200 U (*)
FL5EN0362	Uracil-DNA Glycosylase (1 U/ μ L)	1000 U (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

