

**NucleoSpin gDNA Clean-up XS****Macherey Nagel**

NucleoSpin gDNA Clean-up XS - post clean-up e concentrazione di piccole quantità di grossi frammenti di DNA - raccomandato in particolar modo per piccole quantità di DNA estratto da campioni forensi - basso volume di eluitzione ( $>6 \mu\text{L}$ ): DNA altamente concentrato - metodo più facile e veloce delle unità di filtrazione per microdialisi: 20 min/6 preps - tecnologia: membrana di silice - formato: mini colonna - campione:  $400 \mu\text{L}$  di DNA in soluzione acquosa ( $2 \mu\text{g}$  DNA) - dimensione frammenti: 100 bp- circa 50 kbp - recupero: 60-70% -  $A_{260}/A_{280}$  1,8-1,9 - volume di eluizione 6-10  $\mu\text{L}$  - capacità legante: 3  $\mu\text{g}$  - colonne NucleoSpin® gDNA Clean-up XS, provette per campionamento (2 mL), tamponi.



| Code      | Description                | Packaging |
|-----------|----------------------------|-----------|
| FC140904E | Kit NucleoSpin Clean Up XS |           |
| FC140904L | Kit NucleoSpin Clean Up XS | 50 prep   |
| FC140904N | Kit NucleoSpin Clean Up XS |           |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

