

**NucleoSpin RNA XS Clean-up****Macherey Nagel**

NucleoSpin RNA XS Clean-up - clean-up altamente efficiente di campioni di RNA pre-purificati o da reazioni enzimatiche - completa rimozione degli inibitori della RT-PCR - alto recupero di RNA e concentrazione - lavora su volumi fino a 300  $\mu$ L - volume di eluizione di 5  $\mu$ L - tecnologia: membrana di silice - formato mini colonna, XS design - campione 300  $\mu$ L di soluzioni di RNA contenenti 90  $\mu$ g di RNA - dimensione dei frammenti >200 nt - recupero tipico: 58-95% - A260-A280 1,9-2,1 - tempo di preparazione: 20 min/6 preps - capacità legante 110  $\mu$ g - comprende: colonna NucleoSpin® RNA Clean-up XS con provette, provette di campionamento (2 mL), provette di campionamento (1,5 mL), tamponi.



| Code      | Description                    | Packaging |
|-----------|--------------------------------|-----------|
| FC140903E | Kit NucleoSpin RNA XS Clean-up | 10 PREP.  |
| FC140903L | Kit NucleoSpin RNA XS Clean-up | 50 PREP.  |
| FC140903N | Kit NucleoSpin RNA XS Clean-up | 250 PREP. |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

