

**NucleoSpin DNA FFPE XS****Macherey Nagel**

NucleoSpin DNA FFPE XS - incrementa la qualità del DNA da campioni in formalina o paraffina, non necessita di xilene - formato: mini colonna - tecnologia: membrana di silice - con sovente per paraffina, inodore - tampone per la dissoluzione dei legami del DNA con la formalina - quantità di campione: #7 sezioni (10 µm) di 250 mm² di area totale - volume di eluizione: 5-30 µL - tempo di preparazione: 70 min/6 prep (esclusa lisi) capacità legante: 50 µg - comprende: colonne NucleoSpin® DNA FFPE XS, provette di campionamento (2 mL), sovente per paraffina, tamponi, proteinasi K.



Code	Description	Packaging
FC140980E	Kit NucleoSpin DNA FFPE XS	kit
FC140980L	Kit NucleoSpin DNA FFPE XS	50 prep
FC140980N	Kit NucleoSpin DNA FFPE XS	250 prep.

Accessories

Code	Description	Packaging
FC1740968	PARAFFIN DESSOLVER BUFFER	25 ML



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

