

**NucleoBond Xtra****Macherey Nagel**

NucleoBond Xtra - purificazione di DNA plasmidico ad alta resa- tecnologia: cromatografia a scambio anionico - formato: colonna - dimensione del vettore < 300 kbp - A260/A280 1,8-1,95.

NucleoBond® Xtra Midi: NucleoBond® Xtra Midi Columns with inserted NucleoBond® Xtra Column Filters, tamponi, RNase A

NucleoBond® Xtra Midi Plus: NucleoBond® Xtra Midi Columns with inserted NucleoBond® Xtra Column Filters, NucleoBond® Finalizers, tamponi, RNase A

NucleoBond® Xtra Maxi: NucleoBond® Xtra Maxi Columns with inserted NucleoBond® Xtra Column Filters, tamponi, RNase A

NucleoBond® Xtra Maxi Plus: NucleoBond® Xtra Maxi Columns with inserted NucleoBond® Xtra Column Filters, NucleoBond® Finalizers Large, tamponi, RNase A



Code	Description	Packaging
FC140410E	Kit NucleoBond Xtra Midi	10 prep
FC140410L	Kit NucleoBond Xtra Midi	50 prep
FC140410M	Kit NucleoBond Xtra Midi	100 prep
FC140412E	Kit NucleoBond Xtra Midi Plus	10 prep.
FC140412L	Kit NucleoBond Xtra Midi Plus	50 prep
FC140414E	Kit NucleoBond Xtra Maxi	10 prep
FC140414L	Kit NucleoBond Xtra Maxi	50 prep
FC140414M	Kit NucleoBond Xtra Maxi	100 prep
FC140416E	Kit NucleoBond Xtra Maxi Plus	10 prep
FC140416L	Kit NucleoBond Xtra Maxi Plus	50 prep
FC140410EC	NucleoBond Xtra Midi Columns colonne	

Accessories

Code	Description	Packaging
FC1740316	Buffer di eluizione ELU, per uso con NucleoBond Xtra Midi/Maxi	600 mL
FC1740375	Buffer di lavaggio, per uso con NucleoBond Xtra Midi/Maxi	1000 mL
FC1740415	Nucleo Bond Xtra Combi Rack	1 pz
FC1740417	Nucleo bond xtra buffer set I	1 set
FC1740563	NUCLEOBOND RACK LARGE	1 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

