

**Microprovette Safe-Lock, PP****Eppendorf**

Microprovette Safe-Lock, PP. - Il tappo di chiusura brevettato elimina ogni pericolo di perdita di reagente- Autoclavabili da aperte (a 121°C, per 20 min.)- Apertura e chiusura facile con una mano- Superfici di scrittura opache- Il tappo, a chiusura rapida tramite gancio di sicurezza, evita un'apertura accidentale (es. durante il riscaldamento).- Massima stabilità meccanica per la centrifugazione- Graduazione- Sono disponibili anche microprovette Safe-Lock Eppendorf Biopur®.- Le microprovette da 0,5 ml sono ottimizzate per le tecniche di PCR, quelle da 0,5 ml e 1,5 ml possono essere centrifugate sino a 30000 xg



Code	Description	Packaging
LLG09409331	Microprovette Safe-Lock, 1,5 ml, in	1000 pz.
LLG09409332	Microprovette Safe-Lock, 1,5 ml, giallo	1000 pz.
LLG09409333	Microprovette Safe-Lock, 1,5 ml, rosso	1000 pz.
LLG09409334	Microprovette Safe-Lock, 1,5 ml, blu	1000 pz.
LLG09409336	Microprovette Safe-Lock, 1,5 ml, ambra	1000 pz.
LLG09409358	Microprovette Safe-Lock, 1,5 ml, colori assortiti	1000 pz.
LLG09409320	Microprovette Safe-Lock, 0,5 ml, in	500 pz.
LLG09409352	Microprovette Safe-Lock, 0,5 ml, blu	500 pz.
LLG09409354	Microprovette Safe-Lock, 0,5 ml, giallo	500 pz.
LLG09409355	Microprovette Safe-Lock, 0,5 ml, rosso	500 pz.
LLG09409356	Microprovette Safe-Lock, 0,5 ml, ambra	500 pz.
LLG09409357	Microprovette Safe-Lock, 0,5 ml, colori assortiti	500 pz.
LLG09409335	Microprovette Safe-Lock, 1,5 ml, verde	1000 pz.
LLG09409345	Microprovette Safe-Lock, 2,0 ml, verde	1000 pz.
LLG09409353	Microprovette Safe-Lock, 0,5 ml, verde	500 pz.
LLG09409341	Microprovette Safe-Lock, 2,0 ml, in	1000 pz.
LLG09409342	Microprovette Safe-Lock, 2,0 ml, giallo	1000 pz.
LLG09409343	Microprovette Safe-Lock, 2,0 ml, rosso	1000 pz.
LLG09409344	Microprovette Safe-Lock, 2,0 ml, blu	1000 pz.
LLG09409346	Microprovette Safe-Lock, 2,0 ml, ambra	1000 pz.



Code	Description	Packaging
LLG09409359	Microprovette Safe-Lock, 2,0 ml, colori assortiti	1000 pz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

