

**Vibromulino MM 400****Retsch**

Questo vibromulino è in grado di macinare un gran numero di materiali quali per esempio: pastiglie, minerali, tessuti, terreno e plastiche, raggiungendo particelle finali delle dimensioni di 5 μm . Un breve raffreddamento del campione posto all'interno del contenitore di macinazione è sufficiente per infragilire il campione quando occorre macinare a temperature criogeniche campioni sensibili. Attenzione! Giare e sfere di macinazione da ordinare separatamente.

- Touch display
- Short grinding times and two grinding stations for high sample throughput
- For dry, wet and cryogenic grinding
- Screwable grinding bowls guarantee loss-free grinding
- Cell disruption for DNA/RNA and protein extraction
- Applications in mechanochemistry thanks to grinding time up to 99 h
- Reproducible results due to calibration
- 9 parameter combinations can be stored



Code	Description	Packaging
LLG04688239	Vibromulino MM 400, 100-240 V, 50/60 Hz	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

