

**NUCLEODUR C18 Pyramid, 3 µm****Macherey Nagel****NUCLEODUR® C18 Pyramid, 3 µm**

Stabile in sistemi eluenti 100% acquosi- USP L1

porosità 110Å, dimensioni particelle 1,8 µm, 3 µm, e 5µm ; carbonio 14%
particelle da 7 e 10 µm per separazioni preparative disponibili su richiesta.

selettività di composti polari d'interesse

eccellente deattivazione. Idonea per LC/MS per la bassa capacità legante.

stabilità pH 1-9

Ideale per: analgesici, penicillina, antibiotici, acidi nucleici, vitamine solubili in
acqua, agenti complessi, acidi organici

colonne analitiche EC, dimensioni particelle 3 µm, 14% C

Ø int. 2.0 mm

Lungh. 125 mm



Code	Description	Packaging
MNL76026030	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 3.0 mm, Lungh. 125 mm	1 pz
MNL76026040	EC 125/4 NUCLEODUR C18 Pyramid, 3 µm EC HPLC column EC 125/4 NUCLEODUR C18 Pyramid, 3 µm length: 125 mm, ID: 4 mm	1 pz
MNL76026046	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 4.6 mm, Lungh. 125 mm	1 pz
MNL76026120	EC 150/2 NUCLEODUR C18 Pyramid, 3 µm EC HPLC column EC 150/2 NUCLEODUR C18 Pyramid, 3 µm length: 150 mm, ID: 2 mm	1 pz
MNL76026130	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 3.0 mm, Lungh. 150 mm	1 pz
MNL76026140	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 4.0 mm, Lungh. 150 mm	1 pz
MNL76026146	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 4.6 mm, Lungh. 150 mm	1 pz
MNL76026220	EC 250/3 NUCLEODUR C18 Pyramid, 3 µm EC HPLC column EC 250/2 NUCLEODUR C18 Pyramid, 3 µm length: 250 mm, ID: 2 mm	1 pz
MNL76026230	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 3.0 mm, Lungh. 250 mm	1 pz
MNL76026240	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 4.0 mm, Lungh. 250 mm	1 pz
MNL76026246	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 4.6 mm, Lungh. 250 mm	1 pz
MNL76026320	EC 50/2 NUCLEODUR C18 Pyramid, 3 µm EC HPLC column EC 50/2 NUCLEODUR C18 Pyramid, 3 µm length: 50 mm, ID: 2 mm	1 pz
MNL76026330	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 3.0 mm, Lungh. 50 mm	1 pz
MNL76026340	EC 50/4 NUCLEODUR C18 Pyramid, 3 µm EC HPLC column EC 50/4 NUCLEODUR C18 Pyramid, 3 µm length: 50 mm, ID: 4 mm	1 pz
MNL76026346	NUCLEODUR C18 Pyramid, 3 µm, Ø int. 4.6 mm, Lungh. 50 mm	1 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

