

**VarioPrep colonne preparative NUCLEODUR C18 Pyramid****Macherey Nagel****VarioPrep colonne preparative NUCLEODUR® C18 Pyramid**

Stabile in sistemi eluenti 100% acquosi- USP L1

porosità 110Å, dimensioni particelle 1,8 µm, 3 µm, e 5µm ; carbonio 14%  
particelle da 7 e 10 µm per separazioni preparative dipsonibili su richiesta.

selettività di composti polari d'interesse

eccellente deattivazione. Idonea per LC/MS per la bassa capacità legante.

stabilità pH 1-9

Ideale per: analgesici, penicillina, antibiotici, acidi nucleici, vitamine solubili in

acqua, agenti complessi, acidi organici

dimensioni particelle 5µm, 14% C

Ø int. 10 mm

Lungh. 250 mm



| Code         | Description                                                                                                                                                          | Packaging |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| MNL762271100 | VP 50/10 NUCLEODUR C18 Pyramid, 5 µm VarioPrep HPLC column<br>VP 50/10 NUCLEODUR C18 Pyramid, 5 µm *** with adjustable<br>Endfitting *** length: 50 mm, ID: 10 mm    | 1 pz      |
| MNL762272100 | VarioPrep colonne preparative NUCLEODUR C18 Pyramid, Ø int. 10<br>mm, Lungh. 250 mm                                                                                  | 1 pz      |
| MNL762272210 | VarioPrep colonne preparative NUCLEODUR C18 Pyramid, Ø int. 21<br>mm, Lungh. 250 mm                                                                                  | 1 pz      |
| MNL762273100 | VP 125/10 NUCLEODUR C18 Pyramid, 5 µm VarioPrep HPLC<br>column VP 125/10 NUCLEODUR C18 Pyramid, 5 µm *** with<br>adjustable Endfitting *** length: 125 mm, ID: 10 mm | 1 pz      |
| MNL762273210 | VP 125/21 NUCLEODUR C18 Pyramid, 5 µm VarioPrep HPLC<br>column VP 125/21 NUCLEODUR C18 Pyramid, 5 µm *** with<br>adjustable Endfitting *** length: 125 mm, ID: 21 mm | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
Reagents operates with  
a Certified Quality  
Management System

