

**Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm****Macherey Nagel****Colonne analitiche EC NUCLEODUR® Sphinx RP, 1.8µm**

selettività distintiva basata sulla superficie bifunzionale - USP L1 e USP L11
dimensione pori 110Å, dimensione particelle 1,8 µm, 3 µm e 5 µm; 14 % C
alta densità di legami covalenti silani per picchi scodati per allargare lo sviluppo del metodo
stabilità pH 1 - 10
adatte per LC/MS grazie alle caratteristiche di basso spurgo
alta riproducibilità e consistenete qualità grazie alle strette procedure QC
range applicativo: antibiotici al chinolone, sulfamidici, xantine, sostituti aromatici
dimensioni particelle 1.8µm, 14% C

Ø int. 2.0 mm
Lungh. 30 mm



Code	Description	Packaging
MNL76082120	Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm, Ø int. 2.0 mm, Lungh. 30 mm	1 pz
MNL76082130	Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm, Ø int. 3.0 mm, Lungh. 30 mm	1 pz
MNL76082140	Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm, Ø int. 4.0 mm, Lungh. 30 mm	1 pz
MNL76082146	Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm, Ø int. 4.6 mm, Lungh. 30 mm	1 pz
MNL76082220	Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm, Ø int. 2.0 mm, Lungh. 50 mm	1 pz
MNL76082230	Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm, Ø int. 3.0 mm, Lungh. 50 mm	1 pz
MNL76082240	Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm, Ø int. 4.0 mm, Lungh. 50 mm	1 pz
MNL76082246	Colonne analitiche EC NUCLEODUR Sphinx RP, 1.8µm, Ø int. 4.6 mm, Lungh. 50 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

