

**VarioPrep colonne preparative NUCLEODUR Sphinx RP****Macherey Nagel****VarioPrep colonne preparative NUCLEODUR® Sphinx RP**

selettività distintiva basata sulla superficie bifunzionale - USP L1 e USP L11
dimensione pori 110Å, dimensione particelle 1,8 µm, 3 µm e 5 µm; 14 % C
alta densità di legami covalenti silani per picchi scodati per allargare lo sviluppo del metodo
stabilità pH 1 - 10
adatte per LC/MS grazie alle caratteristiche di basso spurgo
alta riproducibilità e consistenete qualità grazie alle strette procedure QC
range applicativo: antibiotici al chinolone, sulfamidici, xantine, sostituti aromatici
dimensione particelle 5 µm, 14% C

Ø int. 10 mm
Lungh. 50 mm



Code	Description	Packaging
MNL762372100	VarioPrep colonne preparative NUCLEODUR Sphinx RP, Ø int. 10 mm, Lungh. 50 mm	1 pz
MNL762373100	VarioPrep colonne preparative NUCLEODUR Sphinx RP, Ø int. 10 mm, Lungh. 250 mm	1 pz
MNL762373210	VP 250/21 NUCLEODUR SPHINX RP, 5 µm VarioPrep HPLC column VP 250/21 NUCLEODUR SPHINX RP, 5 µm *** with adjustable Endfitting *** length: 250 mm, ID: 21 mm	1 pz
MNL762375100	VP 125/10 NUCLEODUR SPHINX RP, 5 µm VarioPrep HPLC column VP 125/10 NUCLEODUR SPHINX RP, 5 µm *** with adjustable Endfitting *** length: 125 mm, ID: 10 mm	1 pz
MNL762375210	VP 125/21 NUCLEODUR SPHINX RP, 5 µm VarioPrep HPLC column VP 125/21 NUCLEODUR SPHINX RP, 5 µm *** with adjustable Endfitting *** length: 125 mm, ID: 21	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

