

**VarioPrep colonne preparative NUCLEODUR 100-10 C18 ec, 10 µm****Macherey Nagel****VarioPrep colonne preparative NUCLEODUR® 100-10 C18 ec, 10 µm**

Disponibile con octadecyl a media densità (C18 - USP L1) e octyl (C8 - USP L7) modificati, dimensione pori 110Å, dimensione particelle 3 µm e 5 µm; su richiesta 7 µm, 10 µm, 12 µm, 16 µm, 20 µm, 30 µm e 50 µm per separazioni preparative, per analisi giornaliere di routine e rilevamento per HPLC preparativa, stabilità pH 1 - 9. Contenuto carbone 17.5% C per C18, 10,5% C per C8, alta riproducibilità da lotto a lotto, per applicazioni standard di routine in cromatografia in fase inversa.

Fase ottadecil, 17,5%, dimensione particelle 10µm

Ø int. 10 mm

Lungh. 250 mm



Code	Description	Packaging
MNL762010100	VarioPrep colonne preparative NUCLEODUR 100-10 C18 ec, 10 µm, Ø int. 10 mm, Lungh. 250 mm	1 pz
MNL762010210	VP 250/21 NUCLEODUR 100-10 C18 ec VarioPrep HPLC column VP 250/21 NUCLEODUR 100-10 C18 ec *** with adjustable Endfitting *** length: 250 mm, ID: 21 mm	1 pz
MNL762011100	VarioPrep colonne preparative NUCLEODUR 100-10 C18 ec, 10 µm, Ø int. 10 mm, Lungh. 50 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

