



**Viscosimetro routine Cannon-Fenske, calibrato per misure manuali**

**Xylem Analytics**

Viscosimetro routine Cannon-Fenske in accordo con ISO 3105 e ASTM D446.

- Volume di riempimento: 7 ... 10 ml;
- Lungh. totale: ca. 245 mm.



| Code        | Description   | Packaging |
|-------------|---|-----------|
| LLG09268913 | Viscosimetro routine Cannon-Fenske, tipo 513, K=0,015 | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

