



Viscosimetro routine Cannon-Fenske, calibrato per misure manuali

Xylem Analytics

Viscosimetro routine Cannon-Fenske in accordo con ISO 3105 e ASTM D446.

- Volume di riempimento: 7 ... 10 ml;
- Lungh. totale: ca. 245 mm.



Code	Description	Packaging
LLG09268913	Viscosimetro routine Cannon-Fenske, tipo 513, K=0,015	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

