

**pH/Redox meter pH 3110 Set SM PRO, LLG Premium Line**Xylem Analytics Germany
(WTW)

The LLG Premium Line is our best selection of high quality and proven devices from our top suppliers.

The meter pH 3110 Set SM PRO is ready-to-use for portable pH and Redox measurement.

Convincing quality

- With proven, universal pH electrode Sentix 41
- Robust with sturdy housing made of ABS
- Additional protection due to protective cover SM Pro made of silicone
- Waterproof according to IP 67

Excellent functionality and best user comfort

- Ready-to-use in carrying case
- Large, clear LCD segment display, easy to read thanks to high contrast
- Clear keyboard with only 6 keys made of a gap-free silicone mat
- Keyboard can be operated with gloves due to tactile feedback
- Easy to clean

Highest safety

- Automatic AutoRead function for reproducible and reliable readings
- 3-point calibration with technical or DIN/NIST buffers
- Integrated calibration timer for reminding to calibrate regularly

Scope of supply: pH meter pH 3110 ProfiLine, pH electrode Sentix 41, buffer STP 4 / STP 7 (50 ml each), protective cover SM Pro, KCL solution, measuring beaker (50 ml), electrode holder, carrying case



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG04684148 | pH/Redox meter pH 3110 Set SM PRO, LLG Premium Line, incl. pH electrode Sentix 41 and protective cover SM Pro. | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

