



Tetracloroetilene

Synonyms:

Percloroetilene

C2Cl4	Classification	Danger
Molecular weight	165.83 ONU	1897
CAS	127-18-4 Transport hazard class:	6.1
EEC-N:	204-825-9 Packing	III

Tetracloroetilene > RS - Per spettroscopia ottica

Grado Speciale

Clear, colourless liq. appearance	Conform	-
Identification	Conform	-
Refractive index at 20°C	1.5014 - 1.5074	-
Density d20/4	1.618 - 1.628	-
Boiling point	120.3 - 121.8	°C
Water (K.F.)	<= 100	ppm



Residue on evaporation	<= 10	ppm
Assay (GC)	>= 99.9	%
Free acid	<= 0.0003	meq/g
Free alkali	<= 0.00006	meq/g
Fluorescence (quinine)	-	-
at 365 nm	<= 2	ppb
UV Transmittance	-	-
at 300 nm	>= 75	%
at 310 nm	>= 83	%
at 320 nm	>= 88	%
at 340 nm	>= 88	%
at 360 nm	>= 88	%
at 400 nm	>= 97	%

ID

Size

Packaging

P0682716

1 l

Bottiglia di vetro

P0682721

2.5 l

Bottiglia di vetro

**CARLO ERBA Reagents S.r.l.**VIA R. MERENDI, 22
20007 CORNAREDO (MI)

TEL. +39 02 93 99 190

FAX +39 02 93 991 001

www.carloerbareagents.com*CARLO ERBA Reagents
operates with a Certified
Quality Management
System*