



Carta Test Qualitativa

Macherey Nagel

Carta Test Qualitativa. Per test rapidi.



Code	Description	Packaging
MNL90601	Carta reattiva qualitativa per rame Cuprotesmo, 0.05g Cu su superfici	40 test
MNL90609	Carta Test Qualitativa, Tipo Watesmo	1 pz
MNL90612	Carta Test Qualitativa Cyantesmo: Rotolo lungo 5 m e largo 10 mm, Tipo Phosphatesmo MI, fosfatasi alcalina	50 pz
MNL90721	Carta Test Qualitativa per alluminio/zirconio, 10 mg/L Al ³⁺ - 20mg/L Zr ⁴⁺	100 test
MNL90726	Carta Test Qualitativa per ferro, 10 mg/L Fe ²⁺ o Fe ³⁺	100 test
MNL90729	Carta per test rame 20 mg/L Cu ²⁺	200 strip
MNL90730	Carta Test Qualitativa per nichel, 10 mg/L Ni ²⁺	200 test
MNL90732	Carta Test Qualitativa per argento, 20 mg/L Ag ⁺	200 test
MNL90750	Carta reattiva qualitativa per fluoruro: 20 mg/L F ⁻	200 pz
MNL90751	Cartina Indanthrene gialla 20x70 mm	200 strip
MNL90765	Carta Test Qualitativa INDIPRO per proteine, 50 ug BSA (albumina da siero bovino)	60 test