

Technical Specification

PRODUCT : Methanol + 0.1% v/v formic acid RS For LC/MS

CODE : 414860

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
HPLC Gradient	-	-
at 254 nm	mAU	<=10
UV Transmittance	-	-
at 230 nm	%	>=10
at 254 nm	%	>=90
Assay (GC)	%	>=99.5
Test LC-MS TIC (100-2000m/z) ESI(+)	-	-
Sensitive Impurities (reserpine)	ppb	<=100
Metals content	-	-
Na	ppm	<=2
Ca	ppm	<=0.5
Mg	ppm	<=0.5
K	ppm	<=0.5
Raw material used	-	-
Methanol LC-MS (code 414830)	-	Batch number
Formic acid 98-99% (code 405820)	-	Batch number

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in