

Technical Specification

PRODUCT : SULFURIC ACID d = 1.820 RS - For agroalimentary analysis

CODE : 502020

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
Density at 20°C	-	1.815 - 1.825

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of inappropriate use.