



**Cápsulas para evaporación, DURAN<sup>®</sup>**

- DIN 12336
- Con pico y base plana
- Esterilizables en autoclave



Code	Description	Packaging
LLG09000032	Evaporating basin DURAN®,diam.50 mm,hei	1 pz
LLG09000034	Evaporating basin DURAN,diam. 60 mm,hei	1 pz
LLG09000038	Evaporating basin DURAN®,diam.70 mm,hei	1 pz
LLG09000041	Evaporating basin DURAN®,diam.80 mm,hei	1 pz
LLG09000044	Evaporating basin DURAN®,diam.95 mm,hei	1 pz
LLG09000049	Evaporating basin DURAN®,diam.115 mm,he	1 pz
LLG09000054	Evaporating basin DURAN®,diam.140 mm,he	1 pz
LLG09000059	Evaporating basin DURAN®,diam.190 mmwit	1 pz
LLG09000063	Evaporating basin DURAN®,diam.230 mmwit	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

