



Cápsulas para evaporación, DURAN[®]

DURAN Group

Cápsulas para evaporación, DURAN[®]

-DIN 12336

-Con pico y base plana

-Esterilizables en autoclave



Code	Description	Packaging
LLG09000032	Evaporating basin DURAN®,diam.50 mm,hei	1 pz
LLG09000034	Evaporating basin DURAN,diam. 60 mm,hei	1 pz
LLG09000038	Evaporating basin DURAN®,diam.70 mm,hei	1 pz
LLG09000041	Evaporating basin DURAN®,diam.80 mm,hei	1 pz
LLG09000044	Evaporating basin DURAN®,diam.95 mm,hei	1 pz
LLG09000049	Evaporating basin DURAN®,diam.115 mm,he	1 pz
LLG09000054	Evaporating basin DURAN®,diam.140 mm,he	1 pz
LLG09000059	Evaporating basin DURAN®,diam.190 mmwit	1 pz
LLG09000063	Evaporating basin DURAN®,diam.230 mmwit	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

