

**Accesorios para el mezclador MICROTRON[®] MB 550/MB 800**

Los accesorios de mezcla MBA 125 a 1000 ml de vidrio pueden ser desmontados por completo fácilmente.

Los accesorios de mezcla MBA 2000 y MBA 4000 están hechos de acero inoxidable. Sus robustas tapas de seguridad tienen una abertura de alimentación integrada para la dosificación. Unidad de cuchillos con 4 cuchillas de acero inoxidable, fácil de desmontar para facilitar la limpieza.

Para el funcionamiento del mezclador 550 MB MICROTRON es obligatoria una cubierta de seguridad, pedirla por separado.



Code	Description	Packaging
LLG06241880	Mixing attachment MBA 500 500 ml workin	1 pz.
LLG06241881	Safety hood MBH 500 for MBA 500	1 pz.
LLG06257342	Mixing attachment MBA 125 for MB 550125	1 pz
LLG06257343	Safety hood MBH 125for MBA 125	1 pz
LLG06259417	Mixing attachment MBA 4000 4000 ml work	1 pz.
LLG06260922	Mixing attachment MBA 250 250 ml workin	1 pz
LLG06260923	Mixing attachment MBA 1000 1000 ml work	1 pz.
LLG06260924	Mixing attachment MBA 2000 2000 ml work	1 pz.
LLG06260925	Safety hood MBH 250 for MBA 250	1 pz
LLG06260926	Safety hood MBH 1000 for MBA 1000	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

