



Condensador Peltronic

Condensador Peltronic: control electrónico completamente automático de la temperatura en la superficie refrigerada y del ventilador



Code	Description	Packaging
LLG09880852	Emissions capacitor Peltronicpower supp	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

