

**Cucharas de medición, HDPE****VITLAB****Cucharas de medición, HDPE**

Color natural. Forma cónica con borde de relleno afilado. Aptas para alimentos conforme a la normativa CE nr. 10/2011.



Code	Description	Packaging
LLG07008959	Measuring scoop, cap. 15 mlPE-HD, lengt	1 pz.
LLG07008960	Measuring scoop, cap. 25 mlPE-HD, lengt	1 pz.
LLG07008961	Measuring scoop, cap. 65 mlPE-HD, lengt	1 pz.
LLG07008962	Measuring scoop, cap. 110 mlPE-HD, leng	1 pz.
LLG09201025	Measuring scoop, cap. 150 mlPE-HD, leng	1 pz.
LLG09201030	Measuring scoop, cap. 350 mlPE-HD, leng	1 pz.
LLG09201035	Measuring scoop, cap. 750 mlPE-HD, leng	1 pz.
LLG09201040	Measuring scoop, cap. 1250 mlPE-HD, len	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

