

**LLG-Filter tips 2.0 Low Retention, PP****LLG Labware****LLG-Filter tips 2.0 Low Retention, PP**

The new generation of LLG Labware filter tips with low liquid retention. Available in easy-to-open boxes with color coded racks. Compatible with all common pipette brands. You can find the appropriate compatibility list on [www.llg-labware.com](http://www.llg-labware.com) in the Downloads section.

- Improved optical clarity and pipette fit
- Colored coded rack inserts for easy identification
- Graduated
- Sterile
- DNase and RNase free



Code	Description	Packaging
LLG04668784	LLG TIPS LR2 0.1-10 $\mu$ L clear rack 10x96	960 pz
LLG04668785	LLG TIPS LR2 0.1-10 $\mu$ L XLclear rack 10x96	960 pz
LLG04668786	LLG TIPS LR2 1-20 $\mu$ L clear rack 10x96	960 pz
LLG04668787	LLG TIPS LR2 1-100 $\mu$ L clear rack 10x96	10x96 pz
LLG04668788	LLG TIPS LR2 1-200 $\mu$ L clear rack 10x96	10 x 96 pz
LLG04668789	LLG TIPS LR2 100-1000 $\mu$ L clear rack 8x96	8 x 96



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

