



**Cucharas dosificadoras, acero inoxidable**

**Bochem**

**Cucharas dosificadoras, acero inoxidable**

Fabricadas en acero inoxidable Remanit® 4301.



Code	Description	Packaging
LLG06080286	Dispensing scoops, SS, 2 L, length 365mm	1 pz.
LLG09201219	Scoops, 18/8-steel, length 190 mm	1 pz.
LLG09201222	Scoops, 18/8-steel, length 225 mm	1 pz.
LLG09201232	Scoops, 18/8-steel, length 320 mm	1 pz.
LLG09201235	Scoops, 18/8-steel, length 350 mm	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

