



Cucharas dosificadoras, acero inoxidable

Bochem

Cucharas dosificadoras, acero inoxidable

Fabricadas en acero inoxidable Remanit® 4301.



Code	Description	Packaging
LLG06080286	Dispensing scoops, SS, 2 L, length 365mm	1 pz.
LLG09201219	Scoops,18/8-steel,length 190 mm	1 pz.
LLG09201222	Scoops,18/8-steel,length 225 mm	1 pz.
LLG09201232	Scoops,18/8-steel,length 320 mm	1 pz.
LLG09201235	Scoops,18/8-steel,length 350 mm	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

