




Tubes inserts, stainless steel for Shaking water baths Hydro

Lauda Dr.R.Wobser GmbH &
Co.KG

Tubes inserts, stainless steel for Shaking water baths Hydro

Inserts to be placed on the shaking tray.

For tubes diam. 25 ... 29 mm



Code	Description	Packaging
LLG09905997	Tube racks for 12x50ml Flacon tubesfor	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

