

**Tapones de rosca, alta temperatura, PBT****DURAN Group****Tapones de rosca, alta temperatura, PBT**

Rojos. Rosca DIN. Con junta revestida de PTFE. Resistentes a temperaturas de hasta 140°C. Esterilizables en autoclave



Code	Description	Packaging
LLG09209014	Screw cap, GL 14PBT, closed, red,with P	1 pz
LLG09209015	Screw caps, GL 14PBT, closed, red,wi/10	10 pz.
LLG09209018	Screw cap, GL 18PBT, closed, red,with P	1 pz
LLG09209019	Screw caps, GL 18PBT, closed, red,wi/10	10 pz.
LLG09209025	Screw cap, GL 25PBT, closed, red,with P	1 pz
LLG09209026	Screw caps, GL 25PBT, closed, red,wi/10	10 pz.
LLG09209032	Screw caps,PMT,red,with PTFE-coated sea	1 pz
LLG09209033	Screw caps,PMT,red,with PTFE-coated /10	10 pz.
LLG09209044	Screw caps,PMT,red,with PTFE-coated /10	10 pz.
LLG09209045	Screw caps,PMT,red,with PTFE-coated sea	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

