

**Tapones de rosca, alta temperatura, PBT****DURAN Group****Tapones de rosca, alta temperatura, PBT**

Rojos. Rosca DIN. Con junta revestida de PTFE. Resistentes a temperaturas de hasta 140°C. Esterilizables en autoclave



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09209014 | Screw cap, GL 14PBT, closed, red,with P | 1 pz |
| LLG09209015 | Screw caps, GL 14PBT, closed, red,wi/10 | 10 pz. |
| LLG09209018 | Screw cap, GL 18PBT, closed, red,with P | 1 pz |
| LLG09209019 | Screw caps, GL 18PBT, closed, red,wi/10 | 10 pz. |
| LLG09209025 | Screw cap, GL 25PBT, closed, red,with P | 1 pz |
| LLG09209026 | Screw caps, GL 25PBT, closed, red,wi/10 | 10 pz. |
| LLG09209032 | Screw caps,PMT,red,with PTFE-coated sea | 1 pz |
| LLG09209033 | Screw caps,PMT,red,with PTFE-coated /10 | 10 pz. |
| LLG09209044 | Screw caps,PMT,red,with PTFE-coated /10 | 10 pz. |
| LLG09209045 | Screw caps,PMT,red,with PTFE-coated sea | 1 pz |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

