

**Tapones de rosca, alta temperatura, PBT****DURAN Group****Tapones de rosca, alta temperatura, PBT**

Rojos. Rosca DIN. Con junta revestida de PTFE. Resistentes a temperaturas de hasta 140°C.  
Esterilizables en autoclave



Code	Description	Packaging
LLG09209014	Screw cap, GL 14PBT, closed, red,with P	1 pz
LLG09209015	Screw caps, GL 14PBT, closed, red,wi/10	10 pz.
LLG09209018	Screw cap, GL 18PBT, closed, red,with P	1 pz
LLG09209019	Screw caps, GL 18PBT, closed, red,wi/10	10 pz.
LLG09209025	Screw cap, GL 25PBT, closed, red,with P	1 pz
LLG09209026	Screw caps, GL 25PBT, closed, red,wi/10	10 pz.
LLG09209032	Screw caps,PMT,red,with PTFE-coated sea	1 pz
LLG09209033	Screw caps,PMT,red,with PTFE-coated /10	10 pz.
LLG09209044	Screw caps,PMT,red,with PTFE-coated /10	10 pz.
LLG09209045	Screw caps,PMT,red,with PTFE-coated sea	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

