

**Tapones de rosca, PP****DURAN Group****Tapones de rosca, PP**

- Con retén labial
- Resistentes a temperaturas de hasta 140°C
- Sólo aptos para frascos de laboratorio DURAN®



Code	Description	Packaging
LLG09072400	Screw capPP, GL 45, yellow,with lip sea	1 pz
LLG09072401	Screw capPP, GL 45, green,with lip seal	1 pz
LLG09072402	Screw capPP, GL 45, grey,with lip seal	1 pz
LLG09209062	Screw caps,PP,blue,GL 25pack of 10/10	10 pz.
LLG09209063	Screw caps,PP,blue,GL 32pack of 10/10	10 pz.
LLG09209064	Screw caps,PP,blue,GL 45pack of 10/10	10 pz.
LLG09209065	Screw caps,PP,blue,GL 25	1 pz
LLG09209067	Screw caps,PP,blue,GL 32	1 pz
LLG09209069	Screw caps,PP,blue,GL 45	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

