



**Tapones de rosca, PP**

**DURAN Group**

**Tapones de rosca, PP**

- Con retén labial
- Resistentes a temperaturas de hasta 140°C
- Sólo aptos para frascos de laboratorio DURAN®



Code	Description	Packaging
LLG09072400	Screw capPP, GL 45, yellow,with lip sea	1 pz
LLG09072401	Screw capPP, GL 45, green,with lip seal	1 pz
LLG09072402	Screw capPP, GL 45, grey,with lip seal	1 pz
LLG09209062	Screw caps,PP,blue,GL 25pack of 10/10	10 pz.
LLG09209063	Screw caps,PP,blue,GL 32pack of 10/10	10 pz.
LLG09209064	Screw caps,PP,blue,GL 45pack of 10/10	10 pz.
LLG09209065	Screw caps,PP,blue,GL 25	1 pz
LLG09209067	Screw caps,PP,blue,GL 32	1 pz
LLG09209069	Screw caps,PP,blue,GL 45	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

