



Anillos de vertido, PP

DURAN Group

Anillos de vertido, PP

Esterilizables en autoclave hasta 140°C con vapor de agua.



Code	Description	Packaging
LLG09209075	Filling ring PP green, GL45	1 pz
LLG09209076	Filling ring PP grey, GL45	1 pz
LLG09209077	Filling ring PP yellow, GL45	1 pz
LLG09209078	Filling rings PP green GL45 pack 10	10 pz.
LLG09209079	Filling rings PP grey GL45 pack 10	10 pz.
LLG09209080	Filling rings PP yellow GL45 pack 10	10 pz.
LLG09209081	Pouring ring, PP, GL 32blue	1 pz
LLG09209082	Pouring rings PP GL 32 blue pack 10	10 pz.
LLG09209083	Pouring ring, PP, GL 45blue	1 pz
LLG09209084	Pouring rings PP GL 45 blue pack 10	10 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

