



Pinzas con puntas revestidas de silicona, acero inoxidable

Kern & Sohn

Pinzas con puntas revestidas de silicona, acero inoxidable



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG06054809 | Stainless steel tweezer 250 mm with sili | 1 pz. |
| LLG09900790 | Stainless steel tweezers 105 mm with si | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

