



Spicemen Jar, DURAN

DWK Life Sciences

Spicemen Jar, DURAN®

Cylindrical. With fitted stoppers.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG09007108 | Spicemen Jar, DURAN®, 100 ml, Grinding 34.5 mm, Diam. 52 mm, Height 135 mm | 1 pz |
| LLG09007116 | Spicemen Jar, DURAN®, 300 ml, Grinding 45 mm, Diam. 69 mm, Height 163 mm | 1 pz |
| LLG09007121 | Spicemen Jar, DURAN®, 600 ml, Grinding 50 mm, Diam. 81 mm, Height 214 mm | 1 pz |
| LLG09007122 | Spicemen Jar, DURAN®, 750 ml, Grinding 60 mm, Diam. 90 mm, Height 240 mm | 1 pz |
| LLG09007126 | Spicemen Jar, DURAN®, 1200 ml, Grinding 60 mm, Diam. 100 mm, Height 253 mm | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

