



Tapones KIM-KAP para frascos de cultivo, PP

Kimble Kontes

Tapones KIM-KAP™ para frascos de cultivo, PP
Esterilizable en autoclave. Codificado con color.



Code	Description	Packaging
LLG06235364	Closure Kim-KapPP, natual, 16 mm, / 1000	1 pezzo
LLG06235365	Closure Kim-KapPP, natural, 20 mm, /500	1 pezzo
LLG06256899	Closure Kim-KapPP, natural, 13 mm,/ 1000	1 pezzo
LLG06256900	Closure Kim-KapPP, natural, 18 mm,/ 1000	1 pezzo
LLG06256904	Closure Kim-KapPP, red, 16 mm, pac/ 1000	1 pezzo



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

