

**Tubo, PTFE****BOLA****Tubo, PTFE**

Altamente resistente a los productos químicos y a temperaturas de -200°C a 250°C.



| Code        | Description                              | Packaging |
|-------------|--|-----------|
| LLG06291372 | PTFE tubing 0.50 x 1.00mm 0.25 wall, pac | 10 pz.    |
| LLG06291373 | PTFE tubing 0.50 x 1.60mm 0.50 mm wall,  | 10 pz.    |
| LLG06291374 | PTFE tubing 0.80 x 1.60mm 0.40 mm wall,  | 10 pz.    |
| LLG06291375 | PTFE tubing 1.00 x 1.60mm 0.30 mm wall,  | 10 pz.    |
| LLG06291376 | PTFE tubing 1.00 x 2.00mm 0.50 mm wall,  | 10 pz.    |
| LLG06291377 | PTFE tubing 1.60 x 3.20mm 0.80 mm wall,  | 10 pz.    |
| LLG06291378 | PTFE tubing 2.00 x 3.00mm 0.50 mm wall,  | 10 pz.    |
| LLG06291379 | PTFE tubing 2.00 x 4.00mm 1.00 mm wall,  | 10 pz.    |
| LLG06291380 | PTFE tubing 4.00 x 6.00mm 1.00 mm wall,  | 10 pz.    |
| LLG06291381 | PTFE tubing 6.00 x 8.00mm 1.00 mm wall,  | 10 pz.    |
| LLG06291382 | PTFE tubing 8.00 x 10.00mm 1.00 mm wall, | 10 pz.    |
| LLG09205608 | PTFE tubing 0.50 x 1.00mm0.25 wall, p.m  | 1 pz.     |
| LLG09205609 | PTFE tubing 0.50 x 1.60mm0.50 mm wall p  | 1 pz.     |
| LLG09205616 | PTFE tubing 0.80 x 1.60mm0.40 mm wall,   | 1 pz.     |
| LLG09205619 | PTFE tubing 1.00 x 1.60mm0.30 mm wall,   | 1 pz.     |
| LLG09205621 | PTFE tubing 1.00 x 2.00mm0.50 mm wall,   | 1 pz.     |
| LLG09205630 | PTFE tubing 1.50 x 2.10mm0.30 mm wall,   | 1 pz.     |
| LLG09205633 | PTFE tubing 1.50 x 3.50mm1.00 mm wall,   | 1 pz.     |
| LLG09205634 | PTFE tubing 1.60 x 3.20mm0.80 mm wall,   | 1 pz.     |
| LLG09205644 | PTFE tubing 2.00 x 3.00mm0.50 mm wall,   | 1 pz.     |



| Code        | Description                            | Packaging |
|-------------|--|-----------|
| LLG09205645 | PTFE tubing 2.00 x 4.00mm1.00 mm wall, | 1 pz.     |
| LLG09205657 | PTFE tubing 3.00 x 4.00mm0.50 mm wall, | 1 pz.     |
| LLG09205658 | PTFE tubing 3.00 x 5.00mm1.00 mm wall, | 1 pz.     |
| LLG09205663 | PTFE tubing 4.00 x 5.00mm0.50 mm wall, | 1 pz.     |
| LLG09205665 | PTFE tubing 4.00 x 6.00mm1.00 mm wall, | 1 pz.     |
| LLG09205668 | PTFE tubing 5.00 x 6.00mm0.50 mm wall, | 1 pz.     |
| LLG09205672 | PTFE tubing 6.00 x 8.00mm1.00 mm wall, | 1 pz.     |
| LLG09205678 | PTFE tubing 7.00 x 9.00mm1.00 mm wall, | 1 pz.     |
| LLG09205680 | PTFE tubing 8.00 x 10.00mm             | 1 pz.     |
| LLG09205691 | PTFE tubing 10.00 x 12.00mm            | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

