

**Tubo, PTFE****BOLA****Tubo, PTFE**

Altamente resistente a los productos químicos y a temperaturas de -200°C a 250°C.



Code	Description	Packaging
LLG06291372	PTFE tubing 0.50 x 1.00mm 0.25 wall, pac	10 pz.
LLG06291373	PTFE tubing 0.50 x 1.60mm 0.50 mm wall,	10 pz.
LLG06291374	PTFE tubing 0.80 x 1.60mm 0.40 mm wall,	10 pz.
LLG06291375	PTFE tubing 1.00 x 1.60mm 0.30 mm wall,	10 pz.
LLG06291376	PTFE tubing 1.00 x 2.00mm 0.50 mm wall,	10 pz.
LLG06291377	PTFE tubing 1.60 x 3.20mm 0.80 mm wall,	10 pz.
LLG06291378	PTFE tubing 2.00 x 3.00mm 0.50 mm wall,	10 pz.
LLG06291379	PTFE tubing 2.00 x 4.00mm 1.00 mm wall,	10 pz.
LLG06291380	PTFE tubing 4.00 x 6.00mm 1.00 mm wall,	10 pz.
LLG06291381	PTFE tubing 6.00 x 8.00mm 1.00 mm wall,	10 pz.
LLG06291382	PTFE tubing 8.00 x 10.00mm 1.00 mm wall,	10 pz.
LLG09205608	PTFE tubing 0.50 x 1.00mm0.25 wall, p.m	1 pz.
LLG09205609	PTFE tubing 0.50 x 1.60mm0.50 mm wall p	1 pz.
LLG09205616	PTFE tubing 0.80 x 1.60mm0.40 mm wall,	1 pz.
LLG09205619	PTFE tubing 1.00 x 1.60mm0.30 mm wall,	1 pz.
LLG09205621	PTFE tubing 1.00 x 2.00mm0.50 mm wall,	1 pz.
LLG09205630	PTFE tubing 1.50 x 2.10mm0.30 mm wall,	1 pz.
LLG09205633	PTFE tubing 1.50 x 3.50mm1.00 mm wall,	1 pz.
LLG09205634	PTFE tubing 1.60 x 3.20mm0.80 mm wall,	1 pz.
LLG09205644	PTFE tubing 2.00 x 3.00mm0.50 mm wall,	1 pz.



Code	Description	Packaging
LLG09205645	PTFE tubing 2.00 x 4.00mm1.00 mm wall,	1 pz.
LLG09205657	PTFE tubing 3.00 x 4.00mm0.50 mm wall,	1 pz.
LLG09205658	PTFE tubing 3.00 x 5.00mm1.00 mm wall,	1 pz.
LLG09205663	PTFE tubing 4.00 x 5.00mm0.50 mm wall,	1 pz.
LLG09205665	PTFE tubing 4.00 x 6.00mm1.00 mm wall,	1 pz.
LLG09205668	PTFE tubing 5.00 x 6.00mm0.50 mm wall,	1 pz.
LLG09205672	PTFE tubing 6.00 x 8.00mm1.00 mm wall,	1 pz.
LLG09205678	PTFE tubing 7.00 x 9.00mm1.00 mm wall,	1 pz.
LLG09205680	PTFE tubing 8.00 x 10.00mm	1 pz.
LLG09205691	PTFE tubing 10.00 x 12.00mm	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

