



## Agitadores centrífugos, PTFE

### Agitadores centrífugos, PTFE

Con un eje de acero inoxidable recubierto de PTFE.



Code	Description	Packaging
LLG09197176	Stirrer PTFE 350x6 mm, 50 mm Ø	1 pz.
LLG09197177	Centrifugal stirrers, 450 mmPTFE, 90 mm	1 pz.
LLG09197178	Centrifugal stirrers, 350 mm PTFE, 90 m	1 pz.
LLG09197179	Centrifugal agitating wave450x90 mm Ø,	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6ª Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

