

**Tapones de vidrio, DURAN[®]****DURAN Group****Tapones de vidrio, DURAN[®]**

De vidrio DURAN[®]. Tipo de vidrio 1 / vidrio neutro según USP, EP y JP.

Esterilizables en autoclave.



Code	Description	Packaging
LLG09010180	Glass caps,DURAN®,for neck diam. 31 mm	1 pz
LLG09010185	Glass caps,DURAN®,for neck diam. 46 mm	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

