

**Tapones de vidrio, DURAN[®]****DURAN Group****Tapones de vidrio, DURAN[®]**

De vidrio DURAN[®]. Tipo de vidrio 1 / vidrio neutro según USP, EP y JP. Esterilizables en autoclave.



Code	Description	Packaging
LLG09010180	Glass caps,DURAN®,for neck diam. 31 mm	1 pz
LLG09010185	Glass caps,DURAN®,for neck diam. 46 mm	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

