



Matraces para cultivo, vidrio DURAN[®], cuello recto

DURAN Group

Matraces para cultivo, vidrio DURAN[®], cuello recto

De vidrio DURAN[®]. Matraz Erlenmeyer. De cuello recto. Adecuados para tapones metálicos. Tipo de vidrio 1 / vidrio neutro según USP, EP y JP. Con código de retorno (identificación de lote). Certificado disponible via internet. Esterilizables en autoclave.



Code	Description	Packaging
LLG09010472	Culture flasks,DURAN®,Erlenmeyer shape,	1 pz
LLG09010473	Culture flasks,DURAN®,Erlenmeyer shape,	1 pz
LLG09010474	Culture flasks,DURAN®,Erlenmeyer shape,	1 pz
LLG09010475	Culture flasks,DURAN®,Erlenmeyer shape,	1 pz
LLG09010476	Culture flasks,DURAN®,Erlenmeyer shape,	1 pz
LLG09010477	Culture flasks,DURAN®,Erlenmeyer shape,	1 pz
LLG09010478	Culture flasks,DURAN®,Erlenmeyer shape,	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

