



Matraces para cultivo, vidrio DURAN[®], cuello recto

DURAN Group

Matraces para cultivo, vidrio DURAN[®], cuello recto

De vidrio DURAN[®]. Matraz Erlenmeyer. De cuello recto. Adecuados para tapones metálicos. Tipo de vidrio 1 / vidrio neutro según USP, EP y JP. Con código de retorno (identificación de lote). Certificado disponible via internet. Esterilizables en autoclave.



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09010472 | Culture flasks,DURAN®,Erlenmeyer shape, | 1 pz |
| LLG09010473 | Culture flasks,DURAN®,Erlenmeyer shape, | 1 pz |
| LLG09010474 | Culture flasks,DURAN®,Erlenmeyer shape, | 1 pz |
| LLG09010475 | Culture flasks,DURAN®,Erlenmeyer shape, | 1 pz |
| LLG09010476 | Culture flasks,DURAN®,Erlenmeyer shape, | 1 pz |
| LLG09010477 | Culture flasks,DURAN®,Erlenmeyer shape, | 1 pz |
| LLG09010478 | Culture flasks,DURAN®,Erlenmeyer shape, | 1 pz |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

