




## Tapones de goma

Deutsch & Neumann

**Tapones de goma**  
Gris. DIN 12871. Sólidos.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09230203 | Rubber stoppers, 3,5 x 6,5 x 15 mm high | 1 pz.     |
| LLG09230205 | Rubber stoppers, 5 x 9 x 20 mm high     | 1 pz.     |
| LLG09230208 | Rubber stoppers, 8 x 12 x 20 mm high    | 1 pz.     |
| LLG09230210 | Rubber stoppers, 10,5 x 14,5 x 20 mm hi | 1 pz.     |
| LLG09230212 | Rubber stoppers, 12,5 x 16,5 x 20 mm hi | 1 pz.     |
| LLG09230214 | Rubber stoppers, 14 x 18 x 20 mm high   | 1 pz.     |
| LLG09230217 | Rubber stoppers, 17 x 22 x 25 mm high   | 1 pz.     |
| LLG09230218 | Rubber stoppers, 18 x 24 x 30 mm high   | 1 pz.     |
| LLG09230221 | Rubber stoppers, 21 x 27 x 30 mm high   | 1 pz.     |
| LLG09230223 | Rubber stoppers, 23 x 29 x 30 mm high   | 1 pz.     |
| LLG09230226 | Rubber stopper, 26 x 32 x 30 mm high    | 1 pz.     |
| LLG09230229 | Rubber stoppers, 29 x 35 x 30 mm high   | 1 pz.     |
| LLG09230231 | Rubber stoppers, 31 x 38 x 35 mm high   | 1 pz.     |
| LLG09230236 | Rubber stoppers, 36 x 44 x 40 mm high   | 1 pz.     |
| LLG09230241 | Rubber stoppers, 41 x 49 x 40 mm high   | 1 pz.     |
| LLG09230247 | Rubber stopper, 47 x 55 x 40 mm high    | 1 pz.     |
| LLG09230250 | Rubber stoppers, 50,5 x 59,5 x 45 mm hi | 1 pz.     |
| LLG09230256 | Rubber stoppers, 56 x 65 x 45 mm high   | 1 pz.     |
| LLG09230260 | Rubber stoppers, 60 x 70 x 50 mm high   | 1 pz.     |
| LLG09230264 | Rubber stoppers, 64,5 x 75,5 x 55 mm hi | 1 pz.     |



| Code        | Description                            | Packaging |
|-------------|----------------------------------------|-----------|
| LLG09230271 | Rubber stoppers, 71 x 83 x 60 mm high  | 1 pz.     |
| LLG09230279 | Rubber stoppers, 79 x 92 x 65 mm high  | 1 pz.     |
| LLG09230287 | Rubber stoppers, 87 x 100 x 65 mm high | 1 pz.     |
| LLG09230294 | Rubber stoppers, 94 x 107 x 65 mm high | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

