




**Tapones de goma**

**Deutsch & Neumann**

**Tapones de goma**  
Gris. DIN 12871. Sólidos.



Code	Description	Packaging
LLG09230203	Rubber stoppers, 3,5 x 6,5 x 15 mm high	1 pz.
LLG09230205	Rubber stoppers, 5 x 9 x 20 mm high	1 pz.
LLG09230208	Rubber stoppers, 8 x 12 x 20 mm high	1 pz.
LLG09230210	Rubber stoppers, 10,5 x 14,5 x 20 mm hi	1 pz.
LLG09230212	Rubber stoppers, 12,5 x 16,5 x 20 mm hi	1 pz.
LLG09230214	Rubber stoppers, 14 x 18 x 20 mm high	1 pz.
LLG09230217	Rubber stoppers, 17 x 22 x 25 mm high	1 pz.
LLG09230218	Rubber stoppers, 18 x 24 x 30 mm high	1 pz.
LLG09230221	Rubber stoppers, 21 x 27 x 30 mm high	1 pz.
LLG09230223	Rubber stoppers, 23 x 29 x 30 mm high	1 pz.
LLG09230226	Rubber stopper, 26 x 32 x 30 mm high	1 pz.
LLG09230229	Rubber stoppers, 29 x 35 x 30 mm high	1 pz.
LLG09230231	Rubber stoppers, 31 x 38 x 35 mm high	1 pz.
LLG09230236	Rubber stoppers, 36 x 44 x 40 mm high	1 pz.
LLG09230241	Rubber stoppers, 41 x 49 x 40 mm high	1 pz.
LLG09230247	Rubber stopper, 47 x 55 x 40 mm high	1 pz.
LLG09230250	Rubber stoppers, 50,5 x 59,5 x 45 mm hi	1 pz.
LLG09230256	Rubber stoppers, 56 x 65 x 45 mm high	1 pz.
LLG09230260	Rubber stoppers, 60 x 70 x 50 mm high	1 pz.
LLG09230264	Rubber stoppers, 64,5 x 75,5 x 55 mm hi	1 pz.



Code	Description	Packaging
LLG09230271	Rubber stoppers, 71 x 83 x 60 mm high	1 pz.
LLG09230279	Rubber stoppers, 79 x 92 x 65 mm high	1 pz.
LLG09230287	Rubber stoppers, 87 x 100 x 65 mm high	1 pz.
LLG09230294	Rubber stoppers, 94 x 107 x 65 mm high	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

