




Tiras de temperatura irreversibles de 8 pasos testoterm®

Testo

**Tiras de temperatura irreversibles de 8 pasos testoterm®**

Tiras para 8 temperaturas. 5 tipos cubren el rango de medición de 37 a 260 °C, cada tira muestra 8 temperaturas diferentes. Cada tira mide 50 x 18 mm.



Code	Description	Packaging
LLG09246002	Temperature recording strips, 51 x 18/10	10 pz.
LLG09246003	Temperature recording strips, 51 x 1/10	10 pz.
LLG09246004	Temperature recording strips, 51 x 1/10	10 pz.
LLG09246005	Temperature recording strips, 51 x 1/10	10 pz.
LLG09246006	Temperature recording strips, 51 x 1/10	10 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

